
Scientific Instrumentation and Control with Web of Things

Vignesh Vaidyanathan (M.Sc. - Software Developer)

About Myself


- Software developer @ CALA, LMU, Garching, Germany
 - B.Sc. Electrical Engineering (2016) (India)
 - M.Sc. Physics (2023) (TUM)

- Centre for Advanced Laser Applications
 - Petawatt Peak Power Laser used for Ion Acceleration




Macchi, Andrea. "A Review of Laser-Plasma Ion Acceleration." arXiv: Plasma Physics (2017): 2, Fig. 1

Requirement of Data-Acquisition



Implementation - Data acquisition customized to experiment




Ultrafast Lasers & High Power Lasers

Optical Cavity ~100% Highly reflective mirror Flashlamp (pump source) Crystal Titanium-Sapphire Optical resonator

Fig. Optical Cavity


Pulse Characteristics

- Picosecond or femtosecond pulses
 (1 second to approximately 31.69 million years)
- Pulses per second 1Hz, 10Hz, also, order of Millions (80 Million pulses per second for 80MHz)



Atomic Clocks, Spectroscopy

High Power Lasers




More often - Custom built to research or application requirements

- Ion Acceleration/Ionoacoustics
- Fundamental Physics

Coding/Implementation

- WoT Interaction Affordances (Example Oscilloscope):
 - Properties
 - Time Resolution, Time Interval, Trigger Level, Range
 - Actions
 - Start Acquisition, Stop Acquisition
 - Events
 - Measurement Data Points or Measurement "Ready"
- Sometimes a state machine
 - If camera is capturing video, don't allow (GUI) client to change exposure time (property write)
 - If acquisition is already running, don't start acquisition again (action invoke)
- Code Organization/Design
 - Object Oriented



Cloudo, CC0, via Wikimedia Commons

Properties



```
1 from hololinked.server import Thing
 2 from hololinked.server.properties import Integer
 3 from pyueye import ueye
   class UEyeCamera(Thing):
 6
       def set exposure(self, value):
           cdbl_in = ueye.double(value)
 8
           ret = ueye.is_Exposure(self.handle, ueye.IS_EXPOSURE_CMD_SET_EXPOSURE, cdbl_in,
 9
                                   ueye.sizeof(cdbl_in))
10
           assert return_code_OK(self.handle, ret)
11
12
       def get_exposure(self):
13
           cdbl out = ueye.double()
14
15
           ret = ueye.is_Exposure(self.handle, ueye.IS_EXPOSURE_CMD_GET_EXPOSURE, cdbl_out,
16
                                   ueye.sizeof(cdbl_out))
17
           assert return_code_OK(self.handle, ret)
           return cdbl_out.value
18
```



```
1 from hololinked.server import Thing
 2 from hololinked.server.properties import Number
 3 from pyueye import ueye
 5 class UEyeCamera(Thing):
 6
       def set_exposure(self, value):
 8
           . . .
 9
       def get_exposure(self):
10
11
12
                            GET
                                                 PUT
                                                                DELETE
       exposure = Number(fget=get_exposure, fset=set_exposure, fdel=None,
13
               doc="exposure time of the camera",
14
               bounds=(0, None), observable=True, unit='ms', state=["ON"])
15
```

```
cam = UEyeCamera()
cam.exposure = 5000 # ms
```

• • •

```
1 from hololinked.server import Thing
2 from hololinked.server.properties import Number
3
4 class UEyeCamera(Thing):
5
6 exposure_time = Number(
7 bounds=(0.0, None),
8 inclusive_bounds=(False, True)
9 doc="Exposure time for image in milliseconds",
10 fget=get_exposure,
11 fset=set_exposure,
12 unit='ms',
13 observable=True,
14 ) # type: float
15
```

```
1 "properties" : {
      "exposure time": {
         "title": "exposure_time",
         "description": "Exposure time for image in milliseconds",
         "unit": "ms",
         "type": "number",
         "forms": [
             "href": "https://example.com/stretcher-gitter/exposure-time",
10
             "op": "readproperty",
11
             "htv:methodName": "GET",
12
             "contentType": "application/json"
13
           },
15
             "href": "https://example.com/stretcher-gitter/exposure-time",
             "op": "writeproperty",
17
             "htv:methodName": "PUT",
             "contentType": "application/json"
           },
             "href": "https://example.com/stretcher-gitter/exposure-time/change-event",
21
22
             "op": "observeproperty",
23
             "htv:methodName": "GET",
24
             "contentType": "text/plain",
25
             "subprotocol": "sse"
         ],
         "observable": true,
29
         "exclusiveMinimum": 0
```

Actions



```
1 from hololinked.server import Thing, action
 3 class UEyeCamera(Thing):
 4
       @action()
 5
       def start video(self):
           """Start continuous (free-running) acquisition. Use get_next_image() to retrieve images.""
           ret = ueye.is_CaptureVideo(self.handle, ueye.IS_DONT_WAIT)
 8
 9
           assert return_code_OK(self.handle, ret)
10
       @action(input_schema={... {'blocking' : {'type' : 'boolean'}}}) POST
11
12
       def snap(self, blocking = False):
           if blocking:
13
14
               self._snap()
15
           else:
16
               threading.Thread(target=self._snap).start()
17
       # not a remotely visible method
18
       def _snap()
19
20
           self.start video()
21
           image, timestamp = self.get_next_image(return_timestamp=True)
           self._last_jpeg = base64.b64encode(self.cast_image(image, format='jpeg'))
22
           self.stop_video()
23
```

Events




```
1 from hololinked.server import Thing, action, Event
 3 class UEyeCamera(Thing):
       snap_completed_event = Event(friendly_name="snap-completed",
 5
                           doc="raised when snap is finished, use for long exposure time snaps")
 6
       # not a remotely visible method
 8
       def _snap()
 9
10
           self.start video()
           image, timestamp = self.get_next_image(return_timestamp=True)
11
           self. last ipeg = base64.b64encode(self.cast_image(image, format='jpeg'))
12
           self.snap_completed_event.push(True)
13
14
           self.stop_video()
```

State Machine




```
1 from hololinked.server import Thing, action, Event
3 class UEyeCamera(Thing):
 4
      snap_completed_event = Event(friendly_name="snap-completed",
 5
                           doc="raised when snap is finished, use for long exposure time snaps")
 6
      # not a remotely visible method
8
9
      def_snap()
          self.state_machine.set_state('CAPTURE')
10
11
          self.start_video()
          image, timestamp = self.get_next_image(return_timestamp=True)
12
          self._last_jpeg = base64.b64encode(self.cast_image(image, format='jpeg'))
13
14
          self.snap_completed_event.push(True)
          self.stop_video()
15
          self.state_machine.set_state('ON')
16
17
```




Energy Meter Data


ShotTime	Shot	Run	Flags	LocalFlags	WriteTime	Value	MeasurementTimestamp
2024-02-06 20:51:43.349	34	1	NaN	NaN	2024-02-06 20:51:43.419	0.002115	2024-02-06 20:51:43.404
2024-02-06 20:51:43.447	35	1	NaN	NaN	2024-02-06 20:51:43.513	0.002133	2024-02-06 20:51:43.511
2024-02-06 20:51:43.549	36	1	NaN	NaN	2024-02-06 20:51:43.604	0.002188	2024-02-06 20:51:43.602
2024-02-06 20:51:43.648	37	1	NaN	NaN	2024-02-06 20:51:43.725	0.002136	2024-02-06 20:51:43.723
2024-02-06 20:51:43.746	38	1	NaN	NaN	2024-02-06 20:51:43.816	0.002212	2024-02-06 20:51:43.814
2024-02-07 06:37:50.520	342101	1	NaN	NaN	2024-02-07 06:37:50.573	0.002170	2024-02-07 06:37:50.571
2024-02-07 06:37:50.623	342102	1	NaN	NaN	2024-02-07 06:37:50.695	0.002233	2024-02-07 06:37:50.693
2024-02-07 06:37:50.721	342103	1	NaN	NaN	2024-02-07 06:37:50.786	0.002220	2024-02-07 06:37:50.784
2024-02-07 06:37:50.819	342104	1	NaN	NaN	2024-02-07 06:37:50.877	0.002223	2024-02-07 06:37:50.875
2024-02-07 06:37:50.922	342105	1	NaN	NaN	2024-02-07 06:37:50.998	0.002176	2024-02-07 06:37:50.996

How to synchronise GUIs




Promote State Machine State to an observable property





Bind State Machine State to GUI Elements' Props



To auto-apply in React & Svelte

- 1) State Property in context API
- 2) Make components react to change to state (useEffect hook)

Special Thanks

- Prof. Dr. Jörg Schreiber, Associate Professor, LMU
 - Chair for Experimental Physics Medical Physics and Laser Ion Acceleration
- Anna Schmidt, M.Sc.
 - Ph.D. student Ion Acoustics and Laser Ion Acceleration
 - Co-organiser for hardware

Attribution

- Fig. Optical Cavity Lakkasuo, CC BY-SA 3.0 https://creativecommons.org/licenses/by-sa/3.0, via Wikimedia Commons
- Fig. Pulse Train HartmutG, CC BY-SA 3.0 http://creativecommons.org/licenses/by-sa/3.0/, via Wikimedia
- Fig. Picoscope Pico Technology, CC BY-SA 3.0 https://creativecommons.org/licenses/by-sa/3.0, via Wikimedia Commons
- Fig. Beam cross section AdamJones88, CC BY-SA 4.0 https://creativecommons.org/licenses/by-sa/4.0, via Wikimedia Commons https://commons.wikimedia.org/wiki/File:Gaussian Beam Laser Pulsed.png >
- Fig. Light Pulse By Zueignung Own work, CC BY-SA 3.0, https://commons.wikimedia.org/w/index.php?curid=22226803

https://www.flaticon.com/free-icons/backend" Backend icons created by Flat Icons - Flaticon

• <u>https://www.flaticon.com/free-icons/data</u> - Data icons created by Freepik – Flaticon




https://www.flaticon.com/free-icons/control-room - Control room icons created by Freepik - Flaticon

- Fig. Spectrum Analyzer Rohde & Schwarz GmbH & Co. KG, Mühldorfstraße 15, D-81671 München, http://www.rohde-schwarz.com, CC BY 3.0 https://creativecommons.org/licenses/by/3.0, via Wikimedia Commons
- Egekorkan, CC BY-SA 4.0 https://creativecommons.org/licenses/by-sa/4.0, via Wikimedia Commons
- Fig 7. Signal Generator Mu301, CC BY-SA 3.0 https://creativecommons.org/licenses/by-sa/3.0, via Wikimedia Commons
- Fig 8. Arduino Uno By SparkFun Electronics from Boulder, USA Arduino Uno R3, CC BY 2.0, https://commons.wikimedia.org/w/index.php?curid=26785892

Ultrafast Lasers & High Power Lasers

Pulse Characteristics

- Picosecond or femtosecond pulses
 (1 second to approximately 31.69 million years)
- Pulses per second order of Millions
 (80 Million pulses per second for 80MHz)



Atomic Clocks, Spectroscopy

Radio Frequency
Optical Switch (say 80MHz)

Computational "Virtual" Thing

- Properties computed beam diameter, centre of mass (beam pointing)
- Actions start computation, take control action
- Events computation complete

